

Filtration Volumes of Media or Aqueous Solutions from 100 mL to 10 L?

yes

Clarification or Sterile Filtration?

Clarification

Sterilization

Solution difficult to filter and/or volumes > 5 L?

Yes

no

Which inlet connection do you prefer?

Which inlet connection do you prefer?

Which inlet connection do you prefer?

Hose barb

Luer-Lock

Hose barb

Luer-Lock

Hose barb

Luer-Lock

18076-----N

Prefilter & 0.45 µm PES
Hose barb inlet & outlet
Pack of 100 pieces, non-sterile

18072-----D

Prefilter
Luer-Lock inlet
Hose barb outlet
Pack of 10 pieces, non-sterile

18075-----UPN

0.2 µm PES
Luer-Lock inlet
hose barb outlet
Pack of 100 pieces, non-sterile

18091-----D

Prefilter & 0.2 µm PES
Hose barb inlet & outlet
PTFE automatic venting
Filling bell
Pack of 10 pieces, sterile

18068-----D

Prefiler & 0.2 µm PES
Luer-Lock inlet
Hose barb outlet
PTFE automatic venting
Filling bell
Pack of 10 pieces, sterile

18092-----D

Prefiler & 0.2 µm PES
Luer-Lock inlet & outlet
Pack of 10 pieces, sterile

18089-----D

0.2 µm PES
Hose barb inlet & outlet
PTFE automatic venting
Filling bell
Pack of 10 pieces, sterile

18090-----D

0.2 µm PES
Luer-Lock inlet & outlet
Pack of 10 pieces, sterile

18075-----D

0.2 µm PES
Luer-Lock inlet
hose barb outlet
Pack of 10 pieces, sterile